



ThirdLeG Mobility Aids

8/8/7 Manolaya, Seethammal Colony,
1st cross street, Alwarpet Chennai – 18

Phone numbers: **044-42305128, 08056223234, 08754408196**
thirdleg@gmail.com , www.thirdlegmobility.com

“POWER WHEELCHAIR-KKR1”

*Battery Powered,
Joystick Operated Wheelchair in
“Front” and “Rear” Wheel Drive
Versions*

**Front Wheel Drive for use
Predominantly Outdoor and on
Rough Terrains*

** Rear Wheel Drive for use Indoor
and on Paved Surfaces*



Technical Specification

Motor	- 2 x 300 Watt 24 VDC Permanent Magnet
Wheelchair Control	- Multi-Speed – Multi-Direction Joystick Controller
Cruising Current	- 6 Amps – (under normal conditions)
Battery	- 24V (2 x 12 V), 17Ah VRLA Batteries in Series
Frame	- Fully Foldable/ Detachable Frame with detachable Batteries
Maximum speed	- 9.0 Kmph can be decreased to 6.0 Kmph if required.
Brakes	- Electronic Circuit Braking.
Recharging time with Charger	- 2 hours for minimum charge; 12 hours for full charge
Charger Current (Continuous)	- 2.00 Amps DC @ 24VDC (50 Watts)
Replacement of Battery	- 24 months - depending on usage
Mileage per charge	- 3 Hours non-stop or 12.0 Kms - sufficient for a full day
Maximum weight of Rider	- 110Kgs (Maximum permissible weight – 120 Kgs)
Manual operation	- Wheel Lock / unlock for manual propulsion by Attendant