



ThirdLeG Mobility Aids

8/7 Seethammal Colony, First Cross Street, Alwarpet
Chennai-600018,

Contact phone numbers: **044-42305128, 08056223234**
vksthirdleg@gmail.com ; www.thirdlegmobility.com

LIGHTEST AND SMALLEST FOLDING WHEELCHAIR (KALYAN)

Aluminum folded frame (new design),
Quick release aluminum lithium battery
24V*8AH, make the wheelchair lighter,
N.W. (without battery: 16KG). Quantity of
battery can be changed to one piece or two
pieces.

Solar panels, 24V*10W (solar panel is
optional)

Folding size: 570×290×690mm can be put in
a small bag

Maximum speed: 6km/h

Maximum climbing degree: 12°

Motor: 150W*2, and charger time is 5 h



This kind of electric wheelchair
can folded in 3 seconds, and open it also
only in 3 second, and you can folded very
small that you can put it in a bag.

Front castor and anti-tip wheels all have
equipped with shock absorption
function device. Users will feel more
comfortable.

Standard configuration is a piece of lithium
battery, a special lithium battery charger
and a wheelchair frame (with backrest,
footrest, wheels, motor, controller)

Solar energy charge bag and the other lithium battery need to pay for them.

Wheelchair with electromagnetic brake. When users stay on the slope, the wheelchair will not slip.



It is very safety, convenient, comfortable, practical and beautiful.