



# ThirdLeG Mobility Aids

8/7 Seethammal Colony, First Cross Street, Alwarpet  
Chennai-600018,

Contact phone numbers: **044-42305128, 08056223234**  
[vksthirdleg@gmail.com](mailto:vksthirdleg@gmail.com) ; [www.thirdlegmobility.com](http://www.thirdlegmobility.com)

## PATIENT HOIST

### Technical Specification

Overall length of base .....	1080mm (42.0")
Overall width of base (closed).....	645mm (25")
Internal width of base (open).....	840mm (33")
Maximum height of base (clearance required).....	90mm (4.0")
Full boom height.....	2000mm (80")
Height of spreader bar (point of sling attachment)...	450mm-1700mm(18" – 67")
Load capacity of hoist.....	120kg
Load capacity of sling provided.....	120kg
Overall weight of hoist.....	50.0 kg

Motorized Actuator 24VDC or 220VAC

